

AEROSEPTINE 4

READY-FOR-USE LIQUID AIRBORNE DISINFECTANT
FOR THE TREATMENT OF SURFACES

AIRBORNE DISINFECTION



MISTING



■ APPLICATION

■ Health facilities, animal houses, laboratories...

■ PRODUCT BENEFITS

- Allows treatment of all type of floors, surfaces, furniture that cannot be moved
- Food safety
- Mint fragrance

Fb Franklab®
notre expertise l'Ultra-Propreté



READY
FOR USE



AEROSEPTINE 4

READY-FOR-USE LIQUID AIRBORNE DISINFECTANT FOR THE TREATMENT OF SURFACES

INSTRUCTIONS FOR USE

- **Volume to spread:** Depending on the wanted microbial activity (see Table below).
- **Retrieval time:** 2 hours after spreading.
- **Use:**
 - Spread without human presence
 - Apply only on clean and dry surfaces, cleaned up with a compatible and adapted disinfectant-detergent: **DDN Surf (TDS 1190) or FB Spray (TDS 1830)**
 - Stop air conditioning before spreading
 - Ventilate before retrieval

Note: Contact time may be extended from 1 hour to 1 night for practical issues.

DISINFECTING PROPERTIES

Activity	Norms	Spreading	Contact time
Bactericidal	NFT 72-281	1.36cl/m ³	Depending on the size of the room
Yeasticidal	NFT 72-281	1.09cl/m ³	Depending on the size of the room
Virucidal	EN 14476 + A1 <i>Rotavirus</i>	5.43cl/m ³	Depending on the size of the room

PRECAUTIONS OF USE

See (SDS 010455)

FEATURES

- Aspect: **Liquid**
- Pure pH: **7**
- Expiration date: **2 years**

COMPOSITION

- Didecyltrimethylammonium chloride
- Mint fragrance

PACKAGING

- ① Box of 4 x 5L.....ref. **1045520**



①

Hazardous product - Follow precautions of use. Use biocides cautiously. Always read label and product information before use.