

DDN Surf

**NEUTRAL DETERGENT DISINFECTANT FOR
BIOCLEANING FLOORS AND
SURFACES**



■ PRODUCT ADVANTAGES

- Disinfectant detergent with neutral pH
- Active on Clostridium Difficile
- Compatible with a wide range of materials
- Mint fragrance
- Biocide registration: TP 02 and TP 04
- Active on enveloped viruses



SteriFrance

DDN Surf

NEUTRAL DETERGENT DISINFECTANT FOR BIOCLEANING FLOORS AND SURFACES

INSTRUCTIONS FOR USE

■ Preparation of the biocleaning solution :

■ Wear suitable personal protective equipment.

■ For biocleaning floors and surfaces in the medical sector:

Dilute 20 ml of NMS Surf with 8 litres (0.25%) of mains water and homogenise the solution.

■ For biocleaning of floors and surfaces in the food, industrial and contract sectors:

Dilute 80ml of DDN Surf in 8L (1%) of water and homogenize the solution (160ml of

DDN Surf in 8L of water (2%) for a fungicidal activity

■ As part of an epidemic treatment for Clostridium Difficile :

Dilute 40ml of DDN Surf in 4L (1%) of water and homogenize the solution

MICROBIOLOGICAL PROPERTIES

Activities	Standard	Strains	Concentration	Contact time
Bactericidal	EN 13727	<i>Staphylococcus aureus</i> , <i>Pseudomonas aeruginosa</i> , <i>Enterococcus hirae</i>	0,25%	5 min.
	EN 16615		0,25%	15 min.
	EN 1276	<i>Staphylococcus aureus</i> , <i>Pseudomonas aeruginosa</i> , <i>Enterococcus hirae</i> , <i>E. coli</i>	1%	5 min.
	EN 13697		1%	5 min.
	EN 13727	<i>Enterococcus faecium</i> ERV	0,25%	5 min.
Yeasticidal	EN 13624	<i>Candida albicans</i>	0,25%	5 min.
	EN 16615		0,25%	15 min.
	EN 1650		0,5%	5 min.
	EN 13697	1%	5 min.	
	EN 13624	<i>Candida auris</i>	0,25%	5 min.
Fungicidal	EN 13697	<i>Aspergillus Brasillensis</i>	2%	10 min.
	EN 1650		2%	15 min.
Virucidal against enveloped viruses	EN 14476	<i>Vaccinia Virus</i>	0,25%	10 min.
Active on Rotavirus	EN 14476	<i>Rotavirus</i>	0,5%	5 min.
			0,25%	10 min.
Active on PRV	EN 14476	<i>PRV (Virus modèle de l'hépatite B)</i>	0,5%	5 min.
Active on HIV	EN 14476	<i>HIV</i>	0,25%	1 min.
Tuberculocidal	EN 14348	<i>Mycobacterium terrae</i>	1%	30 min.
			0,25%	15 min.
Sporicidal against Clostridium	EN 13704 EN 17126	<i>Clostridium difficile</i>	1%	15 min.

1: The test on virucidal activity against enveloped viruses (Vaccine Viruses) is representative of the enveloped viruses presented in Annex A of standards EN 14476+A2 and EN 16777 and in Annex B of standard EN 17111, such as: Hepatitis B Virus (HBV), Hepatitis C Virus (HCV), Human Immunodeficiency Virus (HIV), Herpesviridae (Herpes virus), Coronavirus, Influenza Virus, Rabies Virus, Rubella Virus, Measles Virus, ...

PRECAUTIONS FOR USE

Refer to the Safety Data Sheet

FEATURES

- Aspect : Liquid
- Color : Green
- pH at 0,25% : 6,6
- pH at 2% : 7,2
- pH at 1% : 7,3 ± 0,50.

COMPOSITION

■ BIOCIDAL ACTIVE SUBSTANCES :

REACTION MASS COMPOSED OF N,N-DIDECYL-N-(2-HYDROXYETHYL)-N-METHYLAMMONIUM PROPIONATE AND N,N-DIDECYL-N-(2-(2-HYDROXYETHOXY)ETHYL)-N-METHYLAMMONIUM PROPIONATE AND N,N-DIDECYL-N-(2-(2-HYDROXYETHOXY)ETHOXY)ETHYL)-N-METHYLAMMONIUM PROPIONATE: 5.4% (W/W).

N-(3-AMINOPROPYL)-N-DODECYLPROPANE-1,3-DIAMINE: 0.45% (W/W) [CAS: 2372-82-9].

PACKAGING

- 1 Box of 250 doses of 20ml.....ref. **1076042**
- 2 Box of 12 x 1Lref. **1076016**
- 3 Box of 4 x 5 Lref. **1076020**



Product for professional use - follow the precautions for use. Use biocides with care. Read the label and product information before use.



FRANKLAB SAS – 40 rue du Chemin Vert 78610 LE PERRY-EN-YVELINES France.
Tél. : + 33 (0)1 39 44 93 40 ■ Fax : + 33 (0)1 39 44 93 41
contact@sterifrance.com ■ www.sterifrance.com

MADE IN FRANCE

FT
1190.L
v18_03.2025