



ENZYMEX L9

ENZYMATIC DISINFECTANT DETERGENT
FOR THE PRE-DISINFECTION OF MEDICAL
DEVICES



■ APPLICATION

- Operating room, endoscopy, ophthalmology and contact lens

■ PRODUCT BENEFITS

- Tri-enzymatic complex (lipase, amylase and protease) combined with effective surfactants
- Instrumentation integrity guaranteed for 72 hours in immersion-bath
- Can be used in ultrasonic vat
- Proven efficiency on biofilm *Pseudomonas aeruginosa*



ENZYMEX L9

ENZYMATIC DISINFECTANT DETERGENT FOR THE PRE-DISINFECTION OF MEDICAL DEVICES

INSTRUCTION FOR USE

■ Preparation of the solution:

- Wear suitable personal protective equipment
- Dilute the appropriate volume of ENZYMEX L9 in 5L water and homogenize the solution (see dilution table)

dose for 5L	Solution concentration		
	0.5%	2%	5%
	25 ml	100 ml	250 ml

■ Use :

- 1- Soiled instruments have to be fully immersed
- 2- Scrub medical devices and use a swab if necessary
- 3- Contact time : according to the antimicrobial activity desired (see disinfecting properties)
- 4- Rinse thoroughly

Renew the bath for each use (or follow the recommendations local mandates in force). According to the French DGS Instruction of December 1, 2011 (excluding endoscopy and ATNC), the bath must be renewed at least every 24 hours and each time it is macroscopically soiled.

MICROBIOLOGICAL PROPERTIES

Activity	Standards	Strains	Concentration	Contact time
Bactericidal	EN 13727	<i>Staphylococcus aureus, Pseudomonas aeruginosa, Enterococcus hirae</i>	0.5%	5 minutes
	EN 14561		0.5%	10 minutes
	EN 13727	<i>Enterococcus faecium VRE</i>	0.5%	10 minutes
Yeasticidal	EN 13624	<i>Candida albicans</i>	0.5%	5 minutes
	EN 14562		0.5%	10 minutes
	EN 13624	<i>Candida auris</i>	0.5%	10 minutes
Virucidal against enveloped viruses (1)	EN 14476	<i>Vaccinia virus</i>	0.5%	10 minutes
	EN 17111			
Active on Herpesvirus	EN 14476	<i>Herpes virus</i>	0.5%	5 minutes
Active on HIV, DBV	EN 14476	<i>DBV virus (Hepatitis C model virus), HIV</i>	0.5%	10 minutes
Tuberculocidal	EN 14348	<i>Mycobacterium terrae</i>	2%	60 minutes
			5%	10 minutes
			2%	60 minutes
	EN 14563	<i>Mycobacterium terrae</i>	5%	10 minutes

(1): The test relating to virucidal activity against enveloped viruses (Vaccine Virus) is representative of the enveloped viruses presented in Annex A of standards EN 14476+A2 and EN 16777 and in Annex B of the standard EN 17111, such as: Hepatitis B virus (HBV), Hepatitis C virus (HCV), Human immunodeficiency virus (HIV), Herpesviridae (Herpes virus), Coronavirus, Influenza virus, rabies, rubella, measles, etc.

PRECAUTION OF USE

Refer to the Safety Data Sheet

FEATURES

- Aspect : Liquid
- Color : Green
- pH at 0.5% : 7.5

COMPOSITION

- Cationic surfactants
- Nonionic surfactants
- Enzymes: lipase, amylase and protease
- Sequestering agents

PACKAGING

- 1 Box of 250 doses of 25 mlref. **1043547**
- 2 Box of 12 x1 L cansref. **1043516**
- 3 Box of 4 x 5 L cans.....ref. **1043520**



Product for professional use - follow the precautions for use. Before use, read the label and product information. Class IIb Medical Device