

FB'Spray

READY-TO-USE PH-NEUTRAL
DISINFECTANT DETERGENT FOR
BIOCLEANING SURFACES



■ PRODUCT ADVANTAGES

- Phosphate-free, EDTA-free, dye-free, fragrance-free
- Bottle with foam dispenser
- Compatible with a wide range of materials thanks to chloride-free formulation
- Biocide registration: TP 02
- Effective against enveloped viruses

FB'Spray

READY-TO-USE PH-NEUTRAL DISINFECTANT DETERGENT FOR BIOCLEANING SURFACES

INSTRUCTIONS FOR USE

■ When using for the 1st time:

- Wear suitable PPE and unlock the gun.

■ Use:

- 1- Spray or evenly distribute FB'Spray directly on the surface to be treated or on a dry non-woven single-use wipe.
- 2- Spread the solution/foam over the surface with the wipe

- 3- Contact time: depending on the antimicrobial activity required
- 4- Do not rinse*.

*Doesn't require rinsing, unless otherwise indicated by the manufacturer of the surface coating treated or for surfaces in contact with foodstuffs.

- 5- Discard the surface
- 6- Repeat if necessary.

MICROBIOLOGICAL PROPERTIES

Activities	Standard	Strains	Contact time
Bactericidal	EN 16615	<i>Staphylococcus aureus</i> , <i>Pseudomonas aeruginosa</i> , <i>Enterococcus hirae</i>	2 min.
	EN 13727		5 min.
	EN 14561		10 min.
	EN 1276	<i>Staphylococcus aureus</i> , <i>Pseudomonas aeruginosa</i> , <i>Enterococcus hirae</i> , <i>E.coli</i>	2 min.
	EN 13697		5 min.
	EN 13727	<i>Enterococcus faecium</i> VRE	
Yeasticidal	EN 16615	<i>Candida albicans</i>	2 min.
	EN 13624		
	EN 1650		5 min.
	EN 13697		10 min.
	EN 14562		5 min.
	EN 13624	<i>Candida auris</i>	
Fungicidal	EN 13624 EN 13697 EN 1650	<i>Aspergillus brasiliensis</i>	15 min.
Virucidal against enveloped viruses ⁽¹⁾		<i>Vaccinia virus</i>	5 min.
Virucidie à spectre limité ⁽²⁾	EN 14476	<i>Norovirus</i> , <i>Adénovirus</i>	15 min.
Actifve on Rotavirus, PRV, BVDV		<i>Rotavirus</i> , <i>PRV virus (model hepatitis B virus)</i> , <i>BVDV virus (model hepatitis C virus)</i>	5 min.
Sporicidie vis à vis de Clostridium difficile	EN 13704	<i>Clostridium difficile</i>	15 min.
	EN 17126		

1: The test for virucidal activity against enveloped viruses (*Vaccine Viruses*) is representative of the enveloped viruses presented in Annex A of standards EN 14476+A2 and EN 16777 and in Annex B of standard EN 17111, such as: Hepatitis B virus (HBV); Hepatitis C virus (HCV); Human immunodeficiency virus (HIV); Herpesviridae (Herpes virus); Coronavirus; Influenza, rabies, rubella, measles viruses; ...

2: The limited-spectrum virucidal activity test is representative of enveloped viruses (Annex A of standards EN 14476+A2 and EN 16777 and Annex B of standard EN 17111), and noroviruses, rotaviruses and adenoviruses.

PACKAGING

- 1 Box of 12 flasks of 500 mL..... réf. **10786160**
- 2 Box of 12 x 750 ml bottles..... réf. **1078640**
- 3 Box of 4 cans of 5 L..... réf. **1078620**

Observe the precautions for use. Use biocides with care. Read the label and product information before use.

PRECAUTIONS FOR USE

Refer to the Safety Data Sheet

FEATURES

- Aspect : Fluid liquid
- pH : 6

COMPOSITION

■ BIOCIDAL ACTIVE SUBSTANCES :

- Reaction mass composed of N,N-didecyl-N-(2-hydroxyethyl)-N-methylammonium propionate and N,N-didecyl-N-(2-(2-hydroxyethoxy)ethyl)-N-methylammonium propionate and N,N-didecyl-N-(2-(2-hydroxyethoxy)ethyl)-N-methylammonium propionate: 0.21% (m/m).

- D-gluconic acid, compound with N,N'-bis(4-chlorophenyl)-3,12-diimino-2,4,11,13-tetraazatetradecanediamidine (2:1): 0.09 % (w/w) [CAS :18472-51-0].

