

# PERALEX 9 HECTO +

PERACETIC ACID-BASED  
READY-TO-USE DISINFECTANT



## ■ APPLICATION

- Operating room, endoscopy, optical equipment

## ■ PRODUCT BENEFITS

- Ready-to-use (without activator)
- High-level disinfection in 10 minutes
- Compatibility with materials used in endoscope disinfection benches, in particular corian
- Formulated without aldehydes: no fixation of proteins

# PERALEX 9 HECTO +

## PERACETIC ACID-BASED READY-TO-USE DISINFECTANT

### INSTRUCTION FOR USE

#### ■ Preparation of the disinfectant solution:

- Wear reusable long-sleeved gloves, goggles, an apron and a mask.
- Pour the required quantity of the ready-to-use solution into a tray so that the equipment to be disinfected is fully immersed

#### ■ Use :

- 1- Immerse the medical devices in the solution, respecting the contact time indicated (for endoscopes, circulate the solution through the channels).
- 2- Immerse medical devices in a rinsing bath using microbiologically controlled water (or sterile water in accordance with local recommendations).
- 3- For endoscopes, circulate the rinsing solution through the channels for at least 5 minutes.
- 4- Wipe with a single-use drape and/or dry with medical air.

Repeat the bath every 14 days maximum. Check the stability of the bath every morning before the first use and then every 4 hours with a PERATEST 500 strip. In France, in accordance with current French regulations, repeat the bath every 7 days maximum. Check the stability of the bath every morning before the first use with a PERATEST 500 test strip.

### MICROBIOLOGICAL PROPERTIES

Activities	Standards	Strains	Contact time
Bactericidal	EN 13727 EN 14561	<i>Staphylococcus aureus</i> , <i>Pseudomonas aeruginosa</i> , <i>Enterococcus hirae</i>	5 minutes
Fungicidal	EN 13624 EN 14562	<i>Candida albicans</i> , <i>Aspergillus brasiliensis</i>	5 minutes
Virucidal	EN 14476 EN 17111	<i>Norovirus</i> , <i>Adenovirus</i> , <i>Poliovirus</i> <i>Norovirus</i> , <i>Adenovirus</i>	5 minutes
Mycobactericidal	EN 14348 EN 14563	<i>Mycobacterium terrae</i> , <i>Mycobacterium avium</i>	10 minutes
Sporicidal against <i>Clostridium Difficile</i>	EN 17126	<i>Clostridium difficile</i>	5 minutes
Sporicidal	EN 17126	<i>Bacillus subtilis</i> , <i>Bacillus cereus</i>	10 minutes

### PRECAUTIONS OF USE

Refer to Safety Data Sheet

### FEATURES

- Aspect : Liquide
- Pur pH : 2,5
- Colour : Colourless

### COMPOSITION

- Peracetic acid
- Hydrogen peroxyde
- Acetic acid

### PACKAGING

- 1 Box of 4 x 5 L.....ref. 1047420



Product for professional use - follow the precautions for use. Before use, read the label and product information. Class IIb medical device