

# Virospray

SPORICIDAL DISINFECTANT  
DETERGENT



## ■ ADVANTAGES AND BENEFITS

- Active from 30 seconds. (EN 16615)
- Wet-wipe cleaning and disinfection of Medical Devices, including endocavitary ultrasound probe surfaces
- High level disinfection
- Cleaning and disinfection of food contact surfaces
- Biocide registration : TP 02/04

# Virospray

## SPORICIDAL DISINFECTANT DETERGENT

### INSTRUCTIONS FOR USE

#### 1<sup>st</sup> use of the spray :

- Wear suitable personal protective equipment and unlock the spray

#### Use :

- Spray VIROSPRAY evenly on a dry non-woven single-use wipe
- Wipe the surface of the medical device to be treated
- Contact time: According to the desired antimicrobial activity

- Do not rinse\*
- Throw away the support
- Repeat the operation if necessary

\*Does not require rinsing unless otherwise specified by the manufacturer of the treated surface coating or for surfaces in contact with skin and mucous membranes.

### MICROBIOLOGICAL PROPERTIES

Activities	Standards	Strains	Contact time		
			Clean conditions	Dirty conditions	
Bactericidal	EN 16615	<i>Staphylococcus aureus</i> , <i>Pseudomonas aeruginosa</i> , <i>Enterococcus hirae</i>	30 seconds	2 min.	
	EN 13727		2 min.		
	EN 14561		5 min.		
	EN 17387	-			
	EN 13697	<i>Staphylococcus aureus</i> , <i>Pseudomonas aeruginosa</i> , <i>Enterococcus hirae</i> , <i>E.coli</i>	-		5 min.
	EN 1276		-		
EN 13727	-				
Yeasticidal	EN 16615	<i>Candida albicans</i>	30 seconds	2 min.	
	EN 13624		5 min.	10 min.	
	EN 14562		5 min.	10 min.	
	EN 17387		-	5 min.	
	EN 1650		-	5 min.	
	EN 13697		-	5 min.	
Fungicidal	EN 13624	<i>Aspergillus brasiliensis</i>	5 min.	10 min.	
	EN 14562		-	5 min.	
	EN 17387		-	5 min.	
Virucidal	EN 14476	<i>Poliovirus</i> , <i>Adenovirus</i> , <i>Norovirus</i>	5 min.	10 min.	
	EN 16777	<i>Rotavirus</i>	-	10 min.	
		<i>Adenovirus</i> , <i>Norovirus</i>	-	5 min.	
Tuberculocidal	EN 14348	<i>Mycobacterium terrae</i>	5 min.	5 min.	
	EN 14563				
Mycobactericidal	EN 14348	<i>Mycobacterium avium</i>	5 min.	10 min.	
Sporicidal against Clostridium Difficile	EN 17126	<i>Clostridium difficile</i>	-	10 min.	
Sporicidal	EN 17126	<i>Bacillus subtilis</i>	10 min.	-	
		<i>Bacillus cereus</i>	15 min.	-	

### PRECAUTIONS OF USE

Refer to Safety Data Sheet

### CHARACTERISTICS

- Aspect : Liquid
- pH : 7

### COMPOSITION

BIOCIDAL ACTIVE SUBSTANCES :

- Ethanol : 42 % (m/m) [CAS : 64-17-5]
- Propan-2-ol : 22 % (m/m) [CAS : 67-63-0]
- N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine : 0.12 % (m/m) [CAS 2372-82-9]

### PACKAGING

- Box of 6 x 1L.....ref. **23140M57**
- Box of 4 x 5L .....ref. **23150M**



Product for professional use - please observe safety precautions. Read the label and product information before use. Use biocides with care. Class IIb medical device.